

PESTICIDE FORMULATION & RESIDUE ANALYTICAL CENTRE, PMD, NIPHM, HYDERABAD

Sr. No. in Scope

NABL / NON NABL

Flow chart for analysis of Cypermethrin in formulation sample

Date of Analysis

S.No.	Step	Execution	Executed By
1.	Sample No.		
2.	Name of Sample		
3.	Procedure		
3.1	Preparation of Internal Standard Solution		
3.1.1	Weigh 2.0 g of Di-cyclohexyl Phthalate (DCP) in 100 ml volumetric flask	g	
3.1.2	Note down the S.No. of balance log book		
3.1.3	Dissolve and dilute up to the mark with toluene		
3.2	Preparation of Standard Solution		
3.2.1	Purity of standard	%	
3.2.2	Weigh 0.375 g a.i. of the standard in a 25 ml volumetric flask	g	
3.2.3	Note down the S.No. of balance log book		
3.2.4	Add 10 ml of internal standard solution (3.1.3)	ml	
3.2.5	Dilute up to the mark with toluene		
3.3	Preparation of Sample Solution		
3.3.1	Note down the percent active ingredient content declared on the sample	%	
3.3.2	Weigh sample so as to contain 0.375 g a.i. in a 25 ml volumetric flask	g	
3.3.3	Note down the S.No. of balance log book		
3.3.4	Add 10 ml of internal standard solution (3.1.3)	ml	
3.3.5	Dilute up to the mark with toluene		
4.	GC Parameters		
4.1	Column: Packed with 3 % Dexil 300 on Chromosorb WHP (100 - 120) mesh		
4.1.1	Length x I.D: 90 cm x 3 mm		
4.2	Gas		
4.2.1	Carrier: Nitrogen: 40 ml/min		
4.2.2	Hydrogen: 45 ml/min		
4.2.3	Air: 450 ml/min		
4.3	Temperature		
4.3.1	Oven: 240°C		
4.3.2	Injector: 270°C		
4.3.3	Detector: 270°C		
4.4	Injection Volume: 1 µl		
5.	Result		
	Sample chromatogram no.		
	Standard chromatogram no.		

Name of the Laboratory : Pesticide Formulation & Residue Analytical Centre, PMD, NIPHM, Hyderabad			
Document No.	:	FC-PF-204	Document Name : Flow chart for analysis of Cypermethrin content, % by mass
Revision No.	:	01	Issue Date : 01/07/2011
Revision Date	:	01/07/2013	Next Revision Date : 01/07/2015
Prepared By		Checked By	Approved By & Issued By
Mrs. T. Sridevi (Deputy Technical Manager)		Mr. C.V. Rao (Technical Manager)	Dr Abhay Ekbote (Director PM & Quality Manager)

PESTICIDE FORMULATION & RESIDUE ANALYTICAL CENTRE, PMD, NIPHM, HYDERABAD

6. CALCULATION:

$$\text{Cypermethrin content, \% by mass} = \frac{A_1 \times A'IS'_2 \times M_1}{A'IS'_1 \times A_2 \times M_2} \times P$$

Where,

- A_1 = Peak area of cypermethrin in the sample solution
 $A'IS'_1$ = Peak area of internal standard in the sample solution
 $A'IS'_2$ = Peak area of internal standard in the standard solution
 A_2 = Peak area of cypermethrin in the standard solution
 M_1 = Mass in 'g' of standard cypermethrin in the standard solution
 M_2 = Mass in 'g' of sample taken for test
 P = Percent purity of cypermethrin standard

Result:

SI.No.	Name of test	Result	Unit	Method of Analysis
1.	Active ingredient		%	12015 - 1987 (Reaffirmed 2002)
Remark / Reference :				

Analyzed by	Name	
	Dated signature	
Checked by	Name	
	Dated signature	

Name of the Laboratory : Pesticide Formulation & Residue Analytical Centre, PMD, NIPHM, Hyderabad					
Document No.	:	FC-PF-204	Document Name	:	Flow chart for analysis of Cypermethrin content, % by mass
Revision No.	:	01	Issue Date	:	01/07/2011
Revision Date	:	01/07/2013	Next Revision Date	:	01/07/2015
Prepared By		Checked By		Approved By & Issued By	
Mrs. T. Sridevi (Deputy Technical Manager)		Mr. C.V. Rao (Technical Manager)		Dr Abhay Ekbote (Director PM & Quality Manager)	