

PESTICIDE FORMULATION & RESIDUE ANALYTICAL CENTRE, PMD, NIPHM, HYDERABAD

Sr. No. in Scope

NABL / NON NABL

Flow chart for analysis of Imidacloprid content in formulation sample

		Date of Analysis	
Sl. No.	Step	Execution	Executed By
1.	Sample No.		
2.	Name of Sample		
3.	Procedure		
3.1	Preparation of Mobile Phase		
3.1.1	Mix acetonitrile and water in the proportion of 80:20 (v/v)		
3.1.2	Pass through membrane filter under vacuum		
3.1.3	Homogenize the mixture using a magnetic stirrer		
3.1.4	Allow to attain room temperature		
3.2	Preparation of Internal Standard Solution		
3.2.1	Weigh 0.03g of Acenaphthene in 100 ml volumetric flask.	g	
3.2.2	<i>Note the serial No. of the balance log book</i>		
3.2.3	Add to it 40 ml of mobile phase (3.1.4)		
3.2.4	Keep it for 30 minute with intermittent shaking		
3.2.5	Dilute up to the mark with mobile phase (3.1.4)		
3.3	Preparation of Standard Solution		
3.3.1	Note the purity of the standard	%	
3.3.2	Weigh 0.1g A.I. of standard in 100 ml volumetric flask	g	
3.3.3	<i>Note the serial No. of the balance log book</i>		
3.3.4	Add to it 40 ml of mobile phase (3.1.4)		
3.3.5	Keep it for 30 minute with intermittent shaking		
3.3.6	Dilute up to the mark with mobile phase (3.1.4)		
3.3.7	Pipette out 10ml of solution (3.3.6) to a 100 ml volumetric flask	ml	
3.3.8	Dilute up to the mark with mobile phase (3.1.4)		
3.3.9	Pipette out 20ml of solution (3.3.8) to a 100 ml volumetric flask	ml	
3.3.10	Add 20ml of internal standard solution (3.2.5)	ml	
3.3.11	Dilute up to the mark with mobile phase (3.1.4)		

Name of the Laboratory :				Pesticide Formulation & Residue Analytical Centre, PMD, NIPHM, Hyderabad			
Document No.	:	FC-PF-215	Document Name	:	Flow chart for analysis of Imidacloprid content, % by mass		
Revision No.	:	01	Issue Date	:	1/07/2011		
Revision Date	:	11/11/2013	Next Revision Date	:	11/11/2015		
Prepared By		Checked By		Approved & Issued By			
Mrs. T. Sridevi (Deputy Technical Manager)		Mr. C.V. Rao (Technical Manager)		Dr. Abhay Ekbote (Director PM & Quality Manager)			

PESTICIDE FORMULATION & RESIDUE ANALYTICAL CENTRE, PMD, NIPHM, HYDERABAD

3.4	Preparation of Sample		
3.4.1	Note the percent active ingredient content declared on sample	%	
3.4.2	Weigh 0.1g A.I. of sample in 100 ml volumetric flask	g	
3.4.3	<i>Note the serial No. of the balance log book</i>		
3.4.4	Add to it 40 ml of mobile phase (3.1.4)		
3.4.5	Keep it for 30 minute with intermittent shaking		
3.4.6	Dilute up to the mark with mobile phase (3.1.4)		
3.4.7	Pipette out 10ml of solution (3.4.6) to a 100 ml volumetric flask	ml	
3.4.8	Dilute up to the mark with mobile phase (3.1.4)		
3.4.9	Pipette out 20ml of solution (3.4.8) to a 100 ml volumetric flask	ml	
3.4.10	Add 20ml of internal standard solution (3.2.5)	ml	
3.4.11	Dilute up to the mark with mobile phase (3.1.4)		
4.	HPLC Parameters		
4.1	Column		
4.1.1	Stainless Steel Packed with Novapak C ₁₈		
4.1.2	Length: 250 mm		
4.1.3	I.D.: 4.6 mm		
4.2	Mobile Phase		
4.2.1	Acetonitrile : Water (80 : 20)		
4.2.2	Flow Rate: 1.0 ml/min		
4.3	Detector: UV		
4.4	Wave Length: 278 nm		
4.5	Injection Volume: 20 µl		
5.	Result		
	Sample chromatogram no.		
	Standard chromatogram no.		

Name of the Laboratory : Pesticide Formulation & Residue Analytical Centre, PMD, NIPHM, Hyderabad			
Document No.	:	FC-PF-215	Document Name : Flow chart for analysis of Imidacloprid content, % by mass
Revision No.	:	01	Issue Date : 1/07/2011
Revision Date	:	11/11/2013	Next Revision Date : 11/11/2015
Prepared By		Checked By	
Mrs. T. Sridevi (Deputy Technical Manager)		Mr. C.V. Rao (Technical Manager)	
Approved & Issued By			
Dr. Abhay Ekbote (Director PM & Quality Manager)			

PESTICIDE FORMULATION & RESIDUE ANALYTICAL CENTRE, PMD, NIPHM, HYDERABAD

6. Calculation:

$$\text{Imidacloprid content, \% by mass} = \frac{A_2 \times A_3 \times M_1}{A_1 \times A_4 \times M_2} \times P$$

Where,

- M_1 = Mass in 'g' of imidacloprid standard
 M_2 = Mass in 'g' of sample taken for test
 A_1 = Peak area of imidacloprid in the standard solution
 A_2 = Peak area of imidacloprid in the sample solution
 A_3 = Peak area of internal standard in the standard solution
 A_4 = Peak area of internal standard in the sample solution
 P = Percent purity of imidacloprid in the standard

Result:

Sl. No.	Name of test	Result	Unit	Method of Analysis
1.	Active ingredient		%	1S 15443 : 2004 (Reaffirmed 2009)
Remark / Reference :				

Analyzed by	Name	
	Dated signature	
Checked by	Name	
	Dated signature	

Name of the Laboratory : Pesticide Formulation & Residue Analytical Centre, PMD, NIPHM, Hyderabad					
Document No.	:	FC-PF-215	Document Name	:	Flow chart for analysis of Imidacloprid content, % by mass
Revision No.	:	01	Issue Date	:	1/07/2011
Revision Date	:	11/11/2013	Next Revision Date	:	11/11/2015
Prepared By		Checked By		Approved & Issued By	
Mrs. T. Sridevi (Deputy Technical Manager)		Mr. C.V. Rao (Technical Manager)		Dr. Abhay Ekbote (Director PM & Quality Manager)	