

**PESTICIDE FORMULATION & RESIDUE ANALYTICAL CENTRE**

PMD, NIPHM, HYDERABAD

Sr. No. in Scope

NABL / NON NABL

**Flow Chart for Analysis of Triazophos Content in formulation sample**

		Date of Analysis	
Sl. No.	Step	Execution	Executed By
1.	Sample No.		
2.	Name of Sample		
	Sample Description		
3.	<b>Procedure</b>		
	<b>3.1 Preparation of Mobile Phase</b>		
3.1.1	Mix Iso-octane and 1,4 Dioxane (saturated with 0.15% water) in the proportion of 90:10 (v/v).		
3.1.2	Pass through membrane filter under vacuum.		
3.1.3	Homogenize the mixture using a magnetic stirrer.		
3.1.4	Allow the mixture to attain room temperature.		
	<b>3.3 Preparation of Standard</b>		
3.3.1	Weigh standard equivalent to 110 mg of active ingredient in 50 ml volumetric flask.	g	
3.3.2	Note the Purity of standard	%	
3.3.3	Add to it 10 ml of Toluene and dissolve the standard.		
3.3.4	Dilute up to the mark with mobile phase (3.1.4).		
3.3.5	Volume of solution (3.1.4) pipette out to a 25 ml volumetric flask.	ml	
3.3.6	Dilute up to the mark with solvent mixture (3.1.4).		
	<b>3.4 Preparation of Sample</b>		
3.4.1	Note the Purity of sample	%	
3.4.2	Weigh accurately sample equivalent to 110 mg of active ingredient in 50 ml volumetric flask.	g	
3.4.3	Add to it 10 ml of Toluene to dissolve the sample		
3.4.4	Dilute up to the mark with solvent mixture (3.1.4)		
3.4.5	Volume of solution (3.4.4) pipette out to a 25 ml volumetric flask	ml	
3.4.7	Dilute up to the mark with solvent mixture (3.1.4)		
3.4.10	Dilute to the mark with solvent mixture (3.1.4)		
4.	<b>HPLC Parameters</b>		
	<b>4.1 Column</b>		
4.1.1	Stainless Steel Packed with Silica particle size 5um		
4.1.2	Length: 250 mm		

Name of the Laboratory : <b>Pesticide Formulation &amp; Residue Analytical Centre, PMD, NIPHM, Hyderabad</b>					
Document No.	:	FC-PF-236	Document Name	:	Flow chart for analysis of Triazophos content
Revision No.	:	00	Issue Date	:	25/03/2013
Revision Date	:	--	Next Revision Date	:	25/03/2015
Prepared By		Checked By		Approved & Issued By	
Ms. T. Sridevi (Deputy Technical Manager)		Mr. C.V. Rao (Technical Manager)		Dr Abhay Ekbote (Quality Manager)	

**PESTICIDE FORMULATION & RESIDUE ANALYTICAL CENTRE****PMD, NIPHM, HYDERABAD**

4.1.3	I.D.:	4.6 mm		
<b>4.2 Mobile Phase</b>				
4.2.1	Iso-octane and 1,4 Dioxane (90 : 10)			
4.2.2	Flow Rate:	2.0 ml/min		
<b>4.3 Detector:</b>			UV	
<b>4.4 Wave Length:</b>			254 nm	
<b>4.5 Injection Volume:</b>			20 µl	
5.	<b>Result</b>			
	Sample Chromatogram no.			
	Standard Chromatogram no.			

**6. Calculation:**

$$\text{Triazophos content, \% by mass} = \frac{A_2 \times M_1}{A_1 \times M_2} \times P$$

**Where,** $M_1$  = Mass in 'g' of Triazophos standard $M_2$  = Mass in 'g' of sample taken for test $A_1$  = Peak area of Triazophos in the standard solution $A_2$  = Peak area of Triazophos in the sample solution

P = Percent purity of Triazophos in the standard

**Result:**

Sl. No.	Name of test	Result	Unit	Method of Analysis
1.	Active ingredient		%	IS 14936 : 2001 (Reaffirmed 2007)

Remark / Reference :

Analyzed by	Name	
	Dated signature	
Checked by	Name	
	Dated signature	

Name of the Laboratory : <b>Pesticide Formulation &amp; Residue Analytical Centre, PMD, NIPHM, Hyderabad</b>					
Document No.	:	FC-PF-236	Document Name	:	Flow chart for analysis of Triazophos content
Revision No.	:	00	Issue Date	:	25/03/2013
Revision Date	:	--	Next Revision Date	:	25/03/2015
Prepared By		Checked By		Approved & Issued By	
Ms. T. Sridevi (Deputy Technical Manager)		Mr. C.V. Rao (Technical Manager)		Dr Abhay Ekbote (Quality Manager)	