

**PESTICIDE FORMULATION & RESIDUE ANALYTICAL CENTRE,
PMD, NIPHM, HYDERABAD**

Flow Chart for Wettability Test in WP Formulation

Sl. No.	Step	Date of Analysis	
		Execution	Executed By
1.	Sample No.		
2.	Name of Sample		
3.	Sample Description		
4.	Procedure		
4.1	Weigh accurately about 2 g of the wettable powder.	g	
4.2	<i>Note down the S.No. of the balance log book.</i>		
4.3	Pour the sample (4.1) uniformly and quickly in the beaker containing 100 ml of the standard hard water.		
4.4	Start the stop-watch simultaneously.		
4.5	Note down the time when the whole material is completely submerged into the water.		

Result:

Sr. No.	Name of test	Result	Method of Analysis
1.	Wettability test		IS : 6940 - 1982 (Reaffirmed 1997)

Analyzed by	Name	
	Dated signature	
Checked by	Name	
	Dated signature	

Name of the Laboratory : Pesticide Formulation & Residue Analytical Centre, PMD, NIPHM, Hyderabad			
Document No.	: FC-PF-283	Document Name	: Flow Chart for Wettability test in WP formulation
Revision No.	: 00	Issue Date	: 20/11/2014
Revision Date	: --	Next Revision Date	: 20/11/2016
Prepared By		Checked By	
Approved By		Issued By	
Mrs. T. Sridevi (Deputy Quality Manager)		Dr. Nirmali Saikia (Technical Manager)	
Dr. Abhay Ekbote (Director PM)		Mr. C.V. Rao (Quality Manager)	